TECHNICAL INFORMATION



* Sildent Bite

Sildent Bite is Product by DIN 13903 standard. It a cartridge type a composed of Base and Catalyst.

• INDICATIONS

- Bite registration
- Impression of the occlusal surface.
- Check the occlusal surface.

BENEFITS

- Extremely fast setting time of 1:00 min
- High Hardness
- Outstanding occlusal details and accurate bite registration
- Easy to trim and gring.

PROPERTIES

Physical & Mechanical Property	Standard	
Mixing Time (sec)	Cartridge Type	Auto-mixing
Mixing Ratio	Base : Catalyst	1:1
Working Time (sec)	Over 30sec	
Color Base	Base : Yellow	Catalyst : White
SettingTime (min)	1min (in the mouth)	
Hardness (HD)	More than 20HD	
Linear Dimensional change (%)	Within 0.2%	